

172
Volumen 33(2)

Volumen 33(2) Marzo - Abril 2016







ORIGINALS

Comparison between non-exhaustive critical velocity estimated by heart rate with exhaustive critical velocity and heart rate variability threshold

Comparative study between symmetrical and asymmetrical sports by static structural analysis in adolescent athletes

Rupture of the distal biceps tendon. Why is increasing its incidence? When it should be repaired?

Injury incidence in Brazilian football referees

DOCUMENT OF CONSENSUS

Consensus on use of infiltration in sport.

Document of Consensus of the Spanish Society of Sports Medicine

REVIEW

Recommendations for adult sport athletes with congenital heart diseases



















ANALIZADOR DE LACTATO

Nueva versión del analizador de "Lactate Scout+" **con bluetooth incorporado**. Para supervisiones médicas y

diagnóstico en medicina del deporte. Alternativa flexible y

económica al análisis de laboratorio. Acreditación CE como

aparato de diagnóstico médico.

:: Volumen de muestra necesario: 0,5 microlitros

:: Tiempo de análisis: 10 segundos

:: Peso total: 85 gramos

:: Memoria: 250 resultados con fecha y hora

Pero además, el nuevo Lactate Scout+ incorpora entre otras las siguientes características diferenciadoras:

:: Chip interno que permite la transmisión de datos vía Bluetooth

:: Nueva pantalla LCD con mejor visibilidad

:: Nuevo Menú con fácil acceso

:: Estanqueidad total de la estructura del analizador